


Continuous-Wave Terahertz Systems

Operation Principle, System Configurations & Accessories


TOPTICA Photonics AG

A horizontal bar with a rainbow gradient, representing the spectral range from 190 nm to 0.1 THz.

All Wavelengths.
190 nm - 0.1 THz


What is terahertz radiation?

confidential



- $1 \text{ THz} \leftrightarrow 33 \text{ cm}^{-1} \leftrightarrow 300 \text{ } \mu\text{m} \leftrightarrow 4.1 \text{ meV}$
- Plastic, paper, cardboard, ... transparent to THz waves → **Imaging**
- Many gases and organic solids show THz “fingerprints” → **Spectroscopy**
- Non-ionizing – no health hazards


Frequency-domain platforms



TopSeller: TeraScan

- Tunable DFB lasers
- Fiber-coupled terahertz antennas
- < 10 MHz resolution
- up to 100 dB dynamic range


Time-domain platforms



TopSellers: TeraFlash pro & smart

- Robust ultrafast fiber lasers
- Fiber-coupled terahertz antennas
- TF pro: 6 THz bandwidth, 95 dB PDR
- TF smart: > 4 THz, 1600 pulse traces/s

Superfast screening platform





TopSeller: TeraSpeed

- Robust ultrafast fiber lasers
- Photoconductive emitter + high-bandwidth Schottky receiver
- 500 kS/s (digital), 100 MS/s (analog)


Frequency-domain vs. time-domain THz

confidential



Frequency-Domain (FD) vs. Time-Domain (TD) THz Spectroscopy

| | FD-THz | TD-THz |
|---|---|---|
| Bandwidth | 0.05 – 2.7 THz, limited by laser | 0.1 – 6 THz |
| Peak dynamic range | ~ 100 dB | ~ 100 dB |
| Frequency resolution | 10 MHz | 10 GHz typ. |
| Acquisition time (complete spectrum) | Minutes to hours, depends on resolution and lock-in time | Milliseconds to 1 min. , depends on pulse trace length and # averages |
| Spectral selectivity | Yes | No |





Frequency-domain terahertz generation & detection

Frequency-domain terahertz generation

confidential

- Two lasers @ adjacent frequencies illuminate photomixer
- Applied bias \rightarrow Photocurrent, modulated at beat frequency
- Surrounding antenna emits THz wave
- Terahertz beam is monochromatic
- Tuning the lasers changes THz wavelength



Frequency-domain terahertz detection

confidential


- Second – unbiased – photomixer serves as THz receiver
- THz wave generates time-varying voltage signal $U(t)$
- Laser beat modulates the photoconductance $G(t)$
- Photocurrent $\propto U(t) \times G(t)$

Proportional to THz electric field
And: depends on phase between $U(t)$ and $G(t)$

Lock-in detection:

- Short integration time \Leftrightarrow high measurement speed
(~ 30 s / spectrum)
- Long integration time \Leftrightarrow lower noise floor


Coherent signal detection



Coherent signal detection (I)


confidential

Michelson interferometer



Changing the arm length or the laser frequency “moves” fringes.

cw THz setup




Fiber-optic splitter creates two beam paths.
RX is reference plane.

Coherent signal detection (II)


confidential

Laser temperature change: 1 mK



- No moving parts → cost-effective
- Detection of THz amplitude and phase
- Highly efficient: Dynamic range ~ 100 dB even though $P_{\text{THz}} \sim 100 \mu\text{W}$

cw THz setup



Coherent signal detection (III)


confidential

- Detected photocurrent:

$$I_{RX}(\nu) \propto E_{THz}(\nu) \cdot \cos(2\pi \Delta L \nu / c)$$

- With:

- I_{RX} = receiver photocurrent
- E_{THz} = THz electric field
- ΔL = path length difference to RX = $|L_1 - L_2|$
- ν = THz frequency
- c = Speed of light
- L_1 begins at beam splitter, continues to TX, includes THz path and terminates at RX
- L_2 begins at beam splitter and terminates at RX



- Amplitude of the terahertz electric field is the envelope of detected photocurrent:


$$I_{RX}(\nu) \propto E_{THz}(\nu) \cdot \cos(2\pi \cdot \Delta L \cdot \nu / c)$$

- Determination of the envelope with TeraScan Software:
Maxima and minima are evaluated and adjacent values averaged.

$$Envelope\left(\frac{1}{2}|\nu_{Max} - \nu_{Min}|\right) = \frac{1}{2}(|E_{THz}(\nu_{Max})| + |E_{THz}(\nu_{Min})|)$$

with:

- $\nu_{Max/Min}$ = THz frequency at photocurrent max. / min.




Fringe spacing

confidential

- Fringe spacing is inversely proportional to terahertz path length:

$$\Delta\nu = c/\Delta L$$




- $\Delta\nu$ Fringe spacing (interval between two phase maxima)
 - c Speed of light
 - ΔL Path length difference to RX
- RX fiber is 30 cm longer than TX fiber → extra optical path ~ 43 cm
 - Consequently, terahertz path length of ~ 43 cm corresponds to $\Delta L = 0$
 - Fringe spacing of 2-3 GHz ⇔ path length between 28-33 cm or 53-58 cm
(values may slightly differ due to production tolerances)



Note: Path length is measured from/to virtual focus point, i.e. includes sections within the photomixers (see “back focus distance” in datasheet).

Dynamic range and noise level

confidential



- Dynamic range (DR):

$$DR(v) = 20 \cdot \log(I_{\text{signal}}/I_{\text{noise}})$$

- I_{signal} – Photocurrent in THz receiver:



- depends on frequency
- independent of lock-in integration time

- I_{noise} - Noise level:

- intrinsic property of the photomixer, depends on optical power on the receiver
- Lock-In-Detection: $I_{\text{noise}} \propto 1/\sqrt{\text{integration time}}$
- RMS measurement without terahertz radiation

Single-megahertz frequency steps

confidential



© Univ. of Auckland, New Zealand

Smallest step size: 1 MHz

- ...corresponds to temperature intervals of 40 μ K
- Step size approaches linewidth of DFB lasers

A. Deninger et al., *J. Infrared Milli Terahz Waves* **36:3** (2015) 269

D. Vogt et al., *J. Infrared Milli Terahz Waves* **40:5** (2019) 524



TeraScan 1550



TeraScan 780



TeraBeam



Tuning Range
Extension




Phase Modulation
Extension




GaAs and
InGaAs
Photomixers

TeraScan 780 and TeraScan 1550

confidential



TeraScan 780




TeraScan 1550

Pre-configured systems

- State-of-the-art GaAs or InGaAs photomixers
- Highest bandwidth: **TeraScan 780**
- Highest dynamic range: **TeraScan 1550**

TeraScan components

confidential





A. Roggenbuck et al., *New J. Phys.* **12** (2010) 43017
D. Stanze et al., *J. Infrared Milli Terahz Waves* **32** (2011) 225
A. Deninger et al., *J. Infrared Milli Terahz Waves* **36** (2015) 269

- 1 DFB diode lasers: $\lambda \sim 0.8 \mu\text{m}$ or $\lambda \sim 1.5 \mu\text{m}$
- 2 GaAs or InGaAs photomixers: up to 100 μW output power, peak dynamic range $\sim 100 \text{ dB}$
- 3 High-precision electronics: computerized frequency control, single-MHz frequency steps possible

DFB pro diode lasers

confidential





Tunable DFB pro lasers

- 2 DFB diodes: $\lambda \sim 0.8 \mu\text{m}$ or $\lambda \sim 1.5 \mu\text{m}$
- Computerized frequency control via precise calibration
- Bandwidth 0 – 2 THz (@ 780 nm)
- Fiber-optic beam combination

DFB pro advantages

- Compact, stable lasers
- Temperature accuracy $< 50 \mu\text{K} \leftrightarrow$ excellent frequency resolution
- Customized frequency ranges u.r.



Photomixer = Fast semiconductor + antenna

- GaAs: $0.8\ \mu\text{m}$ → widest tuning range of DFBs
- InGaAs: $1.5\ \mu\text{m}$ → maximum dynamic range
- Max. power ~ $100\ \mu\text{W}$ (InGaAs @ 100 GHz)
- Peak dynamic range > 90 dB (> 100 dB typ.)

Advantages of TOPTICA's photomixers

- Long-standing partnerships with world-leading manufacturers
- Highest dynamic range of any commercial cw-THz system
- Packaged, fiber-pigtailed modules


Control unit: DLC smart

- Digital control electronics with Ethernet interface
- 75% smaller footprint, 50% less weight, compared to previous version
- Improved frequency accuracy due to real-time temperature measurement
- Control unit “recognizes” laser heads connected
- System can be part of a computer network → convenient remote control




Modular product packages

confidential




Tuning Range Extension

confidential



Combination of three DFB lasers, tuning range 600 GHz/laser



- Wavelength coverage from DC to 2.75 THz
- Excellent agreement between measured spectra and HITRAN data



A. Deninger et al., *J. Infrared Milli Terahz Waves*
36:3 (2015) 269.

“Resolution Booster”: Phase and amplitude information @ <10 MHz resolution

- Standard configuration: Spectral resolution depends on phase “fringes”, typical resolution ~ 1 GHz (black symbols in figure)
- Phase modulation extension measures **envelope** of fringes (yellow trace in fig.)
- Technology: 2 fiber stretchers for path length modulation (up to 3 mm @ 1 kHz)
- Twin fiber concept provides enhanced thermal stability




A. Roggenbuck et al., *J. Opt. Soc. Am. B* **29** (2012) 614

TeraScan specifications

confidential

| | TeraScan 780 | TeraScan 1550 |
|--|---|---|
| “In a nutshell” | Wide scan range | High power |
| Optical wavelength | 780 nm | 1550 nm |
| Semiconductor material | GaAs | InGaAs on InP |
| Emitter/Receiver (Article codes) | 2 x photomixer (#EK-000831 / #EK-000832) | High-speed PD / photomixer (#EK-000724 / #EK-000725) |
| Scan range | 0 – 1.8 THz (typ. > 2 THz), without changing lasers | 0 – 1.2 THz Up to 2.7 THz with 3 lasers |
| THz power | Typ. 2 μ W @ 100 GHz, 0.3 μ W @ 500 GHz | 100 μW @ 100 GHz, 10 μW @ 10 GHz |
| Dynamic range @ 300 ms integr. time | > 80 dB @ 100 GHz > 70 dB (typ. 80 dB) @ 500 GHz Typ. 70 dB @ 1000 GHz | > 90 dB (typ. 100 dB) @ 100 GHz > 70 dB (typ. 80 dB) @ 500 GHz Typ. 70 dB @ 1000 GHz |
| Polarization | Circular, ~ 30 dB PER | Linear, ~ 10 dB PER |
| Frequency accuracy | < 2 GHz | |
| Minimum step size | < 10 MHz | |




Optomechanics




Schottky Diodes

Accessories



#BG-002653: Compact rail setup



#OE-000888: Reflection head

- Compact + robust setups
- Rail with 2 mirrors for transmission measurements
- Reflection head with 4 mirrors → beam focus at sample location





#BG-001481: Two mirrors, collimated beam




#BG-001784: Four mirrors, additional focus


- Flexible solutions for transmission measurements
- Precise 3-axis-stages for THz antennas
- Two mirrors: Collimated beam
- Four mirrors: Additional beam focus

| Optomechanics: Specifications | | | | |
|-------------------------------|------------------------------|--------------------------------|--------------------------------|---------------------------------|
| | #BG-002653 (compact rail) | #BG-001481 (2-mirror setup) | #BG-001784 (4-mirror setup) | #OE-000888 (Reflection head) |
| User mode | Transmission | Transmission | Transmission | Reflection |
| Collimating mirrors | Ø 1", focal length 2" | Ø 2", focal length 3" | Ø 2", focal length 3" | Ø 1", focal length 2" |
| Focussing mirrors | -- | -- | Ø 2", focal length 2" | Ø 1", focal length 4" |
| Focus size | -- | -- | ~ 2 mm | ~ 3 mm |





“High-responsivity” model #EK-000933:
Broadband spectrum with GaAs emitter



“High-bandwidth” model #EK-000961:
Detection of individual THz pulses

Schottky receiver: Incoherent detection scheme

- Output signal proportional to THz power, no phase relation with transmitter
- Suitable for spectroscopy, ideal for THz imaging
- Two versions with different amplifier bandwidth available

F. Rettich et al., *J. Infrared Milli THz Waves* **36:7** (2015) 607

M. Yahyapour et al., *IEEE Trans. THz Science Technol.* **6:5** (2016) 670

S. Brinkmann et al., *J. Infrared Milli THz Waves* **38** (2017) 339

Schottky diodes specifications

confidential

| Schottky diodes: Specifications | | |
|---------------------------------|---|--|
| | #EK-000933 ("High Responsivity") | #EK-000961 ("High Bandwidth") |
| Concept | Zero-bias Schottky diode | |
| Terahertz bandwidth | 50 – 1500 GHz | |
| Noise-equivalent power | 7 pW/sqrt(Hz) @ 100 GHz 100 pW/sqrt(Hz) @ 1000 GHz | 70 pW/sqrt(Hz) @ 100 GHz 1000 pW/sqrt(Hz) @ 1000 GHz |
| Typ. responsivity | 22000 V/W @ 100 GHz 1100 V/W @ 1000 GHz | 230 V/W @ 100 GHz 17 V/W @ 1000 GHz |
| Amplifier bandwidth | 10 Hz – 1 MHz | 10 MHz – 4 GHz |
| Applications | Imaging, Spectroscopy | Imaging, THz-based communication, ultrafast dynamic processes |

