

Production and Quality Control Datasheet TeraBeam / TeraScan



1.0 General

Order Number VB.21.02659
Institute / Company Oklahoma State University
Production Date 08.02.2022

2.0 Components Components

System	TeraScan 1550	SN 01107
Laser Head #1	DFB pro BFY THz #LD-1550-0040-DFB-8, 1550 nm	SN 33639
Laser Head #2	DFB pro BFY THz #LD-1550-0040-DFB-8, 1545 nm	SN 33632
Laser Head #3	DFB pro BFY THz #LD-1550-0040-DFB-8, 1530 nm	SN 33633
Control Electronics	DLC smart	SN 015163 Firmware 3.1.1
Emitter	PCA-FD-1550-100-TX-1	SN TC211009
Receiver	PCA-FD-1550-130-RX-1	SN RC191001
Accessories	Schottky receiver Schottky-FD-RX-1 Optomechanics THz Optics / Tr 2M	SN 43-41

3.0 Acceptance Criteria

		Specification	Test Result
Difference Frequency Tuning		2700 GHz	0 .. 3250 GHz
Dynamic Range of Terahertz Power	100 GHz	90 dB	104 dB
@ 300 ms integration time	500 GHz	70 dB	83 dB

Production and Quality Control Datasheet TeraBeam / TeraScan



4.0 Lasers

Frequency Range: 0– 1390 GHz

Settings	THz frequency at startup: 700 GHz		
Laser Connections		DLC smart	fiber port
	Laser SN 33639:	L 1 (Channel 1)	I1
	Laser SN 33632:	L 2 (Channel 2)	I2
Optical Power	Fiber port	O1	O2
	Laser SN 33639:	16.6 mW	16.4 mW
	Laser SN 33632:	16.4 mW	16.6 mW
	Total optical power:	33.0 mW	33.0 mW

Frequency Range: 1130 – 2500 GHz


Settings	THz frequency at startup: 1800 GHz		
Laser Connections		DLC smart	fiber port
	Laser SN 33633:	L 1 (Channel 1)	I1
	Laser SN 33632:	L 2 (Channel 2)	I2
Optical Power	Fiber port	O1	O2
	Laser SN 33633:	16.6 mW	16.4 mW
	Laser SN 33632:	16.4 mW	16.6 mW
	Total optical power:	33.0 mW	33.0 mW

Frequency Range: 1850 – 3250 GHz

Settings	THz frequency at startup: 2600 GHz		
Laser Connections		DLC smart	fiber port
	Laser SN 33633:	L 1 (Channel 1)	I1
	Laser SN 33639:	L 2 (Channel 2)	I2
Optical Power	Fiber port	O1	O2
	Laser SN 33633:	16.6 mW	16.4 mW
	Laser SN 33639:	16.4 mW	16.6 mW
	Total optical power:	33.0 mW	33.0 mW

CAUTION!	Carefully clean the fiber facets before connecting fibers. An inspection under an optical microscope is recommended.
-----------------	---

5.0 Photoconductive Antennas

TX Module	<p>Emitter PCA-FD-1550-100-TX-1</p> <p>InGaAs photodiode with bow-tie antenna</p> <p>Typical terahertz power: 100 μW @ 100 GHz, 10 μW @ 500 GHz</p> <p>Bandwidth: ~ 3 THz</p>	SN TC211009
RX Module	<p>Receiver PCA-FD-1550-130-RX-1</p> <p>InGaAs photomixer with bow-tie antenna</p> <p>Bandwidth: ~ 3 THz</p>	SN RC191001
Package	<p>Housing: cylindrical, \varnothing 25 mm</p> <p>Lens: hyperhemisphere, \varnothing 10 mm, height 6 mm</p> <p>Electrical connector: BNC</p> <p>Fiber: SM/PM fiber, FC/APC connector</p> <p>Fiber length: 1.0 m (emitter), 1.3 m (receiver)</p>	
Beam Properties	<p>Linear THz polarization</p> <div data-bbox="582 931 1302 1232" data-label="Diagram">  </div>	
Operation	<p>Maximum optical power: 35 mW on the chip</p> <p>Maximum bias voltage:</p> <ul style="list-style-type: none"> • Emitter: -2 .. +0.5 V • Receiver: -0.5 .. +0.5 V (only for testing!) 	

CAUTION!	Photoconductive antennas are electrostatic-discharge sensitive devices!
CAUTION!	We recommend to always wear a high-impedance grounding strap for handling.
	Carefully clean the fiber facets before connecting fibers. An inspection under an optical microscope is recommended.

Production and Quality Control Datasheet TeraBeam / TeraScan



Emitter Operation


Optical power: 33 mW @ fiber port O1

Max. DC bias: -2 V DC

Emitter Characteristics

DC photocurrent: -14.3 mA @ -2.0 V

Emitter DC Photocurrent



Receiver Operation

Optical power: 33 mW @ fiber port O2

Max. DC bias: -0.5 V DC (only for testing)


Receiver Characteristics

DC photocurrent: -0.216 mA @ -0.5 V

Noise:

- 119 pA @ 3 ms lock-in integration time
- 9.6 pA @ 300 ms lock-in integration time

Receiver DC Photocurrent



Production and Quality Control Datasheet TeraBeam / TeraScan



Lock-in Detection Settings

Emitter bias: $-0.4 \text{ V DC} \pm 0.9 \text{ V AC}$ (offset \pm modulation amplitude)
 Modulation frequency: 39.7 kHz
 Lock-in phase: 95.3 deg
 Integration time: 3 ms / 300 ms


NOTE!

Set the DC bias offset before adding the AC modulation amplitude.

Terahertz Characteristics

Frequency	100 GHz	500 GHz	1000 GHz
Receiver photocurrent	1464 nA	124 nA	14 nA
Dynamic range @ 300 ms	104 dB	83 dB	63 dB

Spectrum of Terahertz Power



Spectrum acquired with customer optomechanics THz Optics / Tr 2M
 (2 mirrors, focal length = 3")

6.0 Additional Components

Transimpedance Amplifier

PDA-S, gain switch setting: 3

Optical Isolators

Integrated isolator (> 30 dB)
 Additional fiber-inline isolator (50 dB)

Fiber Array


2x2 fiber array, split ratio: 50:50 (#OK-001485)

Production and Quality Control Datasheet TeraBeam / TeraScan




7.0 Instructions and Comments

Sketch of Optical Connections



Sketch of Electrical Connections



8.0 Quality Control

Production Cansu Arpacioğlu Date 08.02.2022

Final Check _____ Date 08.02.2022